# opSenS

**XTender™**

**Catheter**

**Instructions for Use**

**STERILE | EO**

Contents sterile unless enclosed package  
has been opened or damaged.  
Sterilized by Ethylene Oxide (EO).


|                                |          |
|--------------------------------|----------|
| <b>English (Anglais)</b> ..... | <b>2</b> |
| <b>French (Français)</b> ..... | <b>6</b> |

## 1. DEVICE DESCRIPTION

The OpSens XTender catheter is an over-the-wire support/guiding catheter with an atraumatic tip. The catheter is used in combination with guiding catheters or sheaths in interventional procedures to assist with delivery of other interventional devices (including guidewires) and procedural fluids to the coronary and peripheral vasculature. The catheter is offered in four (4) sizes compatible with 5.5Fr, 6Fr, 7Fr and 8Fr guiding catheters or 5Fr, 5.5Fr, 6Fr and 7Fr sheaths. The catheter is used in conjunction with standard length guidewires up to 0.035" (0.89 mm) diameter while accessing the peripheral vasculature and with guidewires up to 0.014" (0.36 mm) diameter while accessing the coronary vasculature. The catheter is 150 cm long comprised of a 125 cm long manipulation shaft connected to a 25 cm distal tube segment. The distal tube segment has a proprietary hydrophilic coating.

The XTender catheter has two (2) radiopaque markers that aid in positioning the catheter. The distal marker is located 2 mm from the distal tip. The proximal marker is located 4 mm from the catheter proximal entry.

The XTender catheter is delivered over a guidewire and through a guiding catheter or sheath resulting in a passageway with an inner diameter that is smaller than that of the guiding catheter or sheath. The passageway inner diameter (or *Catheter I.D.*) is outlined in Table A. Table A also outlines the XTender catheter distal shaft outer diameter (or *Catheter O.D.*). Furthermore, the XTender catheter has a strain relief printed with both the *Catheter I.D.* and *Catheter O.D.*

**CAUTION:** Only physicians trained in vascular surgery, interventional radiology or cardiology, and who have completed training or have experience with intravascular catheters and associated devices should consider using this device.

**CAUTION:** Read the entire Instructions for Use manual prior to using the device.

| Table A: XTender Catheter Specifications |            |                                   |                                   |                     |                     |
|--|------------|-----------------------------------|-----------------------------------|---------------------|---------------------|
| No Article                               | Ref Number | Compatible Guiding Catheter       | Compatible Sheath                 | Catheter I.D.       | Catheter O.D.       |
| 8271-001                                 | F4101      | ≥ 5.5F<br>[ID ≥ 0.066" / 1.68 mm] | ≥ 5F<br>[ID ≥ 0.066" / 1.68 mm]   | 0.052" /<br>1.32 mm | 0.064" /<br>1.63 mm |
| 8271-002                                 | F4102      | ≥ 6F<br>[ID ≥ 0.071" / 1.80 mm]   | ≥ 5.5F<br>[ID ≥ 0.072" / 1.83 mm] | 0.057" /<br>1.45 mm | 0.068" /<br>1.73 mm |
| 8271-003                                 | F4103      | ≥ 7F<br>[ID ≥ 0.078" / 1.98 mm]   | ≥ 6F<br>[ID ≥ 0.079" / 2.01 mm]   | 0.063" /<br>1.60 mm | 0.076" /<br>1.91 mm |
| 8271-004                                 | F4104      | ≥ 8F<br>[ID ≥ 0.088" / 2.24 mm]   | ≥ 7F<br>[ID ≥ 0.092" / 2.34 mm]   | 0.072" /<br>1.83 mm | 0.086" /<br>2.18 mm |

## 2. INDICATIONS FOR USE

The XTender catheter is intended for use in conjunction with guiding catheters or sheaths during coronary and peripheral interventional procedures to guide and support interventional devices, including guidewires, traverse discrete portions of the vasculature, allow for interventional device exchanges and provide a conduit for infusion of saline solution, diagnostic contrast agents and therapeutic agents.

## 3. CONTRAINDICATIONS

This device should not be used in patients who cannot be systemically anticoagulated.

## 4. WARNINGS

- The XTender catheter is supplied STERILE and for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilizing could increase the risk of patient infection and of compromised device performance.
- The catheter should only be manipulated while observing under fluoroscopy.
- If resistance is encountered at any time during the insertion procedure, do not force passage or torque the catheter. Resistance may cause damage to device or vessel. Carefully withdraw the catheter.
- The catheter should only be advanced or withdrawn over a guidewire.
- Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the device or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.
- This device has not been evaluated for use in the neurovasculature.

## 5. PRECAUTIONS

- Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment. Do not use if the package or catheter is damaged since the sterility or integrity of the device may be compromised, increasing the risk of patient infection and device malfunction.
- To avoid damage to hydrophilic coating, do not wipe the catheter surface with dry gauze.
- Do not attempt to pass the catheter through a smaller sized guiding catheter or sheath than indicated in Table A. Damage to the device may occur.
- Only use interventional devices indicated for passage through the XTender catheter I.D. indicated in Table A. Consult the literature provided by the manufacturer of the individual interventional devices for passageway size requirements and usage guidelines.
- To avoid kinking and damaging the catheter, advance slowly in small increments.
- To avoid damage to the catheter or vessel, do not advance or withdraw the device without a guidewire in place.
- Catheter should not be advanced into a vessel having a diameter smaller than the catheter outer diameter (Refer to Table A for XTender catheter O.D.). Damage to the device or vessel may occur. If pressure dampens after inserting the catheter into a vessel, withdraw the catheter until the pressure returns to normal.
- Never advance the XTender catheter more than 15 cm beyond the tip of the guiding catheter or sheath as the catheter may become lodged in the guiding catheter or sheath making it difficult to remove.
- When delivering fluid through the guiding catheter or sheath via the hemostasis valve Y connector, use a hand held syringe.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.
- If strong resistance is encountered during manipulation of the devices, do not force passage. Determine the cause of the resistance before proceeding. If the cause cannot be removed, withdraw all devices simultaneously.
- Use systemic anticoagulation to prevent or reduce clotting, as preclinical testing with this device showed the potential for clot formation in the absence of anticoagulation.

## 6. POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

Complications may occur with the use of any access or support catheter or during any catheterization procedure.

Possible complications associated with this type of procedure include but are not limited to the following:

- Acute or total occlusion of the vessel
- Acute myocardial infarction
- Arterial dissection, perforation, rupture or injury
- Arterial spasm
- Arterial thrombosis
- Catheter fracture with separation and distal embolization
- Death
- Hemorrhage or hematoma

- Hypo/hypertension
- Local or systemic infection including sepsis
- Stroke, air embolism and embolization or fragmentation of thrombotic or atherosclerotic material
- Surgical intervention

## 7. PACKAGING, STERILIZATION AND STORAGE

The device has been sterilized using ethylene oxide (EO) and is supplied sterile and non-pyrogenic. The package label indicates the **Use By** date.

Do not use the device after the **Use By** date. The device sterility and integrity may be compromised and possibly result in patient infection and device malfunction.

The device should be stored in cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures (greater than 60°C or 140°F). Rotate inventory so that device is used prior to the **Use By** date on the package label.

## 8. RECOMMENDED ITEMS

**Each XTender catheter package includes the following:**

- Single-use sterile disposable catheter contained in a coiled storage tube.

**Materials required but not provided are:**

- Guiding catheter or sheath with an inner diameter large enough to accommodate the specific model of XTender catheter intended for use. Refer to Table A for compatibility and catheter dimension information.
- Y-adaptor with hemostasis valve (Tuohy-Borst type).
- Standard length guidewire up to 0.035" (0.89 mm) diameter for accessing the peripheral vasculature and/or guidewire up to 0.014" (0.36 mm) diameter for accessing the coronary vasculature.
- Sterile syringe (with 20 G maximum diameter blunt needle) for flushing catheter inner diameter.
- Sterile heparinized saline for catheter preparation.
- Additional XTender catheters.

## 9. DEVICE PREPARATION

- Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment. Do not use catheter if either the catheter or packaging is damaged or compromised.
- If the packaging is free of damage, carefully open the pouch and introduce the sterile hoop (with catheter) to the sterile field using sterile techniques
- Carefully disengage the catheter handle from the hoop connector and withdraw the catheter from the hoop. Carefully inspect the catheter for any signs of damage. If damaged, please discard and use another catheter.
- The device is coated with a lubricious hydrophilic coating. Prior to inserting the catheter, activate the coating by immersing the catheter in heparinized saline for approximately 30 seconds, or wiping down the catheter shaft with a heparinized saline saturated gauze or sponge. Do not wipe the catheter surface with dry gauze.
- Thoroughly flush the distal tube lumen of the XTender catheter from the distal tip with heparinized saline solution using syringe/needle.

## 10. DEVICE USAGE

- Appropriate anticoagulation and vasodilation therapy should be used.
- Using standard techniques, place the guiding catheter/sheath and guidewire. Ensure a Y-adaptor with hemostasis valve is connected to the guiding catheter or sheath.
- Place the prepared XTender catheter over the previously inserted guidewire by threading the end of the guidewire through the tip of the catheter. Advance the catheter over the guidewire and through the hemostasis valve. Refer to Table A for guidance regarding appropriate catheter selection along with guiding catheter and sheath sizing.
- The catheter can be introduced percutaneously through the guiding catheter or sheath. Do not attempt to pass the device through a smaller sized guiding catheter or sheath introducer than indicated in Table A.
- Using fluoroscopic imaging, carefully advance the catheter to the desired location in the vasculature. If resistance is encountered at any time during the insertion procedure, do not force passage or torque the device. The radiopaque marker on the distal end of the catheter will provide visual fluoroscopic guidance.

The XTender catheter may be advanced up to 15 cm beyond the tip of the guiding catheter or sheath.

- Using fluoroscopic imaging manipulate the guidewire and catheter in order to achieve the desired position. Do not torque the catheter. Also, do not advance or withdraw the catheter without the guidewire in place and extending beyond the distal end of the catheter.

- Perform the catheterization procedure. Follow the individual manufacturer’s instructions for any interventional devices used during the procedure. Consult Table A regarding passageway requirements.
- When removing the catheter, use fluoroscopic imaging and make sure to maintain guidewire position. Remove the catheter prior to removing the guiding catheter or sheath.

**DEVICE EXCHANGE**

When exchanging interventional devices, including guidewires, maintain the XTender catheter position and carefully withdraw the interventional device under fluoroscopic imaging. Without moving the catheter, insert new interventional device through the hemostasis valve, guiding catheter (or sheath) and XTender catheter under fluoroscopic imaging. Do not advance the XTender catheter without the guidewire in place and extending beyond the distal end of the catheter.

Common with all interventional procedures, exchanging/introducing guidewires may result in the guidewire(s) wrapping around each other or an in-dwelling interventional device, including the XTender catheter. This may impede the advancement of interventional devices along the guidewire(s). Use standard interventional techniques to identify and maintain proper guidewire position.

**FLUID INFUSION**

To use the XTender catheter for fluid infusion, maintain the catheter position and withdraw any interventional device other than a guidewire. A guidewire may be left in place during fluid infusion. Ensure the Y-adaptor hemostasis valve is securely tightened around the proximal shaft of the XTender catheter (and possibly the guidewire) to avoid any fluid loss. Attach a hand-held syringe (filled with desired fluid) to the Y connector on the hemostasis valve. If required, a stopcock may be desired for additional fluid control. Inject the desired fluid.

Fluid injections with a hand-held syringe are recommended.

**11. DISPOSAL**

After use, this product is a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable hospital, local, state and federal laws and regulations.

**WARRANTY DISCLAIMER AND LIMITATION OF LIABILITY**

Although the OpSens XTender catheter hereafter referred to as “Product” has been manufactured and designed under carefully controlled conditions, OpSens Inc. has no control over the conditions under which this Product is used.

OpSens Inc., therefore disclaims all warranties, both express and implied, written or oral, with respect to the Product including, but not limited to, any implied warranty as to the condition, quality, durability, performance, merchantability or fitness for a particular purpose.

OpSens Inc. shall not be liable to any person or entity for any medical expenses, any loss or injury to a party’s profits or goodwill or any direct, incidental, consequential, special, punitive or exemplary damages caused by, arising from or related to any use, defect, failure or malfunction of the Product, whether a claim for such damages is based upon warranty, contract, tort or otherwise.

No person has any authority to bind OpSens Inc. to any representation or warranty with respect to the Product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

|   |   |
|---|---|
| <p><b>Manufacturer:</b><br/> <b>OpSens, LLC</b><br/>         750 Boulevard du Parc Technologique<br/>         Québec, QC G1P 4S3<br/> <br/>         Tel: 1-418-781-0333<br/> <br/> <a href="http://www.opsensmedical.com/">www.opsensmedical.com/</a></p> | <p><b>Australian Sponsor:</b><br/> <b>Emergo Australia</b><br/>         201 Sussex Street, Darling Park<br/>         Tower II, Level 20<br/>         Sydney NSW 2000<br/>         Australia<br/>         Tel: +61 2 9006 1662</p> |
|---|---|

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